

VA Application Note No. V-95

Title: Quinine in bitter lemon

Summary: Quinine can be determined by polarography at the DME using Britton-Robinson buffer at pH = 7.0 as supporting electrolyte.

Sample: Bitter lemon (soft drink)

Sample preparation: none

Determination of quinine

Electrolyte Britton-Robinson buffer pH = 7:
 $c(\text{H}_3\text{PO}_4) = 0.04 \text{ mol/L}$
 $+ c(\text{CH}_3\text{COOH}) = 0.04 \text{ mol/L}$
 $+ c(\text{H}_3\text{BO}_3) = 0.04 \text{ mol/L}$,
 adjust pH to 7.0 with NaOH

Measuring solution 10 mL buffer
 + 1 mL sample

Working electrode (WE) MME (Multi-Mode Electrode) 6.1246.020

Auxiliary electrode (AE) Pt 6.0343.000

Reference electrode (RE) Ag/AgCl/KCl (3 mol/L): 6.0728.020 + 6.1245.010


Parameters

Working electrode	DME
Stirrer speed	2000 rpm
Mode	DP
Purge time	300 s
Equilibration time	5 s
Pulse amplitude	50 mV
Start potential	-800 mV
End potential	-1250 mV
Voltage step	6 mV
Voltage step time	0.4 s
Sweep rate	15 mV/s
Peak potential quinine	-1130 mV

Results:	Quinine
	33.5 mg/L

Determination of quinine

**Quinine in Bitter lemon
sample**



Quinin
c = 33.473 mg/l
+/- 0.419 mg/l (1.25%)

