803 Ti Stand

Manual

8.803.8001ML / 2019-09-26

Metrohm AG

CH-9100 Herisau

Switzerland

Phone +41 71 353 85 85

Fax +41 71 353 89 01

info@metrohm.com

www.metrohm.com

803 Ti Stand

Manual

Technical Communication
Metrohm AG
CH-9100 Herisau
techcom@metrohm.com

This documentation is protected by copyright. All rights reserved.

This documentation has been prepared with great care. However, errors can never be entirely ruled out. Please send comments regarding possible errors to the address above.

Table of contents

1	Introduction	1
1.1	Instrument description	1
1.2	About the documentation	1
1.2.1	Symbols and conventions	1
1.3	Safety instructions	2
1.3.1	General notes on safety	2
1.3.2	Electrical safety	2
1.3.3	Working with liquids	3
1.3.4	Flammable solvents and chemicals	3
1.3.5	Recycling and disposal	4
2	Overview of the instrument	5
3	Installation	7
3.1	Setting up the instrument	7
3.1.1	Packaging	7
3.1.2	Checks	7
3.1.3	Location	7
3.2	Mounting the 803 Ti Stand	7
3.3	Connecting the 803 Ti Stand	10
3.4	Mounting the accessories	12
4	Troubleshooting	15
4.1	Problems	15
5	Technical specifications	16
5.1	Rotational speed	16
5.2	Pump	16
5.3	Power supply	16
5.4	Ambient temperature	16
5.5	Dimensions	17
5.6	Material of housing	17
6	Accessories	18

Table of figures

Figure 1	Front 803 Ti Stand	5
Figure 2	Rear 803 Ti Stand	6
Figure 3	Connecting the 803 Ti Stand	10
Figure 4	Connecting the Dosino to the 803 Ti Stand	11
Figure 5	803 Ti Stand with volumetric titration vessel and equipment for auto- matic reagent exchange.	12

There are no parts inside the housing which can be serviced or replaced by the user.

WARNING

Do not fail to pull the power cord out of the power socket before you set up or disconnect electrical plug connections at the rear of the instrument.

CAUTION

WARNING

- Set up the instrument in a well-ventilated location (e.g. fume cupboard).
- Keep all sources of flame far from the workplace.
- Clean up spilled liquids and solids immediately.
- Follow the safety instructions of the chemical manufacturer.

1.3.5 Recycling and disposal

This product is covered by European Directive 2012/19/EU, WEEE – Waste Electrical and Electronic Equipment.

The correct disposal of your old instrument will help to prevent negative effects on the environment and public health.

More details about the disposal of your old instrument can be obtained from your local authorities, from waste disposal companies or from your local dealer.

2 Overview of the instrument


Figure 1 Front 803 Ti Stand

- | | |
|---|---|
| 1 Bore hole
For support rod (6.2016.070), diameter 10 mm, length 40 mm | 2 Housing
With built-in pump and magnetic stirrer |
| 3 LED function display
Lights up, when the stirrer is switched on. | 4 Key
Pressing the key aspirates air out of the aspiration bottle. The vacuum in the suction bottle suctions the liquid out of the KF titration cell and into the suction bottle. |
| 5 Key
Pressing the key pumps air into the solvent bottle. The overpressure in the solvent bottle pushes solvent into the KF titration cell. | |


Figure 2 Rear 803 Ti Stand

1 Connection nipple for PVC tubing
For aspirating the titration vessel content.


2 Connection nipple for PVC tubing
For aspirating solvent.

3 MSB connector
For connecting additional devices (e.g. Dosino)


4 MSB connection cable
For connecting to an analysis device.


- 2** Attach the opening on the bottom of the KF titration stand onto the hexagon screw in the stand plate.

- Slide the 6.2016.070 support rod from above into the opening of the KF titration stand intended for this purpose up to the hexagon screw and tighten.


Mounting the KF titration cell

Install the KF titration cell as follows:

- 1 Screw the 6.2013.010 clamping ring tightly to the support rod.
- 2 Mount the 6.1414.030 vessel lid of the KF titration cell to the support rod. Keep the locking lever pressed down until it can be released at the desired position.
- 3 Fasten the 6.1415.220 (or 6.1415.250) titration vessel with a 6.1903.020 (or 6.1903.030) stirring bar inside on the vessel lid. Fold back the holding bracket upwards while doing so. The markings on the vessel lid and on the plastic ring must be aligned above one another.

- 4** Press the holding bracket downwards in order to fix the titration vessel. The levers of the holding bracket must enclose the pins of the plastic ring on the titration vessel in order to ensure a secure hold.
- 5** Adjust the height of the KF titration cell by pressing the locking lever. It should almost touch the surface of the stirrer. The position can now be fixed by readjusting the clamping ring.

Once the height of the KF titration cell has been adjusted correctly, the entire cell can be raised and swiveled to the right as required by pressing the locking lever.

3.3 Connecting the 803 Ti Stand

Connect the titration stand as follows:

- 1 Exit the control software.
- 2 Connect the connection cable of the titration stand to one of the sockets marked with **MSB** on the rear of the control device.


Figure 3 Connecting the 803 Ti Stand

- 3** Connecting an additional device (e.g. Dosimat or Dosino) to the MSB connector of the titration stand (optional).


Figure 4 Connecting the Dosino to the 803 Ti Stand

- 4** Start the control software.

22 Molecular sieve 6.2811.000

Mounting the equipment for aspirating

- 1 Equip the 6.1602.105 bottle attachment with the 6.1446.040 threaded stopper, the 6.1808.050 tubing adapter and the 6.1446.090 stopper.
- 2 Plug the 6.1623.000 overflow protection from below into the 6.1602.105 bottle attachment (the opening, the 6.1808.050 tubing adapter is plugged in).
- 3 Screw the 6.1602.105 bottle attachment onto the 6.1608.030 waste bottle.
- 4 Cut the 6.1801.120 PVC tubing into two halves and use one to aspirate.
- 5 Fasten the one end of the aspiration tubing to the tubing adapter, the other to the connection nipple for the waste on the rear of the titration stand.
- 6 Insert the 6.1819.050 cannula through the 6.1446.090 stopper into the waste bottle.
- 7 Screw the 6.1805.200 PTFE tubing into the 6.1446.090 stopper.
- 8 Screw the 6.1543.120 aspiration tip to the other end of the PTFE tubing.
- 9 Insert the aspiration tip down to the bottom of the titration vessel to be emptied and fasten it to the titration vessel lid.

NOTICE

The waste bottle should periodically be emptied.

Mounting the equipment for aspirating solvent

Mount the solvent bottle as follows:

- 1** Fill the 6.1609.010 adsorber tube and the 6.2811.000 molecular sieve with cotton.
- 2** Use the other half of the 6.1801.120 PVC tubing previously cut in two pieces for aspirating the solvent. Fasten the one end of the tubing to the lower end of the adsorber tube, the other end to the connection nipple for aspirating the solvent.
- 3** Equip the 6.1602.105 bottle attachment with the 6.1446.040 threaded stopper and the 6.1609.010 adsorber tube.
- 4** Secure the adsorber tube with the 6.2023.020 SGJ clip.
- 5** Insert the 6.1819.030 cannula into the last open hole (with medium-sized diameter) of the 6.1602.105 bottle attachment.
- 6** Screw the 6.1602.105 bottle attachment onto the solvent bottle.
Instead of the 6.1608.023 amber glass bottle, other reagent bottles with GL45 thread can be used (e.g. RIEDEL DE HAEN (1 liter), BAKER). For other bottles another bottle attachment or an additional thread adapter have to be used.
MERCK bottles: 6.1602.110 bottle attachment
FLUKA bottles, RIEDEL DE HAEN (500 ml): 6.1602.100 bottle attachment with 6.1618.000 thread adapter.
- 7** Screw the 6.1543.110 transfer tip to the free end of the 6.1805.200 PTFE tubing, insert into the titration vessel and fasten to the titration vessel lid.

4 Troubleshooting

4.1 Problems

Problem	Cause	Remedy
Solvent is added at the same time during aspiration.	<i>The drying tube is blocked.</i>	▪ Enlarge the hole of the drying tube cover to 2 mm at least. ▪ Refill the drying tube, not packing too tightly, eventually removing some of the molecular sieve.
The pump conveys no liquid.	<i>The add/aspirate key is not pressed tightly enough.</i>	▪ It is usually the case that the bottle attachments are not screwed tightly enough to the reagent and waste bottles. ▪ Check all connections for leak-tightness. ▪ Check whether the 803 Ti Stand has been set up correctly.
The solvent flows on into the titration vessel after the addition, without the add key being pressed.	<i>The solvent bottle is not set up correctly.</i>	Set the solvent bottle up in such a way that the liquid level in the bottle is lower than that in the titration vessel.

<i>Width</i>	106 mm
<i>Height (without stand)</i>	101 mm
<i>Height (with stand)</i>	412 mm
<i>Depth</i>	220 mm
<i>Weight</i>	1100 g

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