



Application Note AN-NIR-073

# Determination of water activity in tablets with the OMNIS NIR Analyzer

Fast, non-destructive measurements performed in seconds

Water activity ( $a_w$ ), the partial vapor pressure of water in a substance divided by the vapor pressure of water in standard state, is expressed either as 0–100% equilibrium relative humidity (ERH) or scaled to 0–1  $a_w$ . It is used to assess the safety, quality, and strength of non-sterile drug pharma products. In compounded preparations,  $a_w$  refers to water that is freely available to participate in reactions (e.g., hydrolysis) or provides an environment that supports microbiological growth. Typically, solid dosage

pharmaceuticals are  $>0.70 a_w$ , indicating that microbial growth is unlikely [1]. Elevated  $a_w$  in powders affects flow, caking, compaction, and strength properties of solid dosage forms and it is used in the study of shelf-life, aging, and packaging requirements. Measuring  $a_w$  in the pharma environment is described in USP<1112> and USP<922> [2]. Dedicated instruments that measure  $a_w$  require up to 30 minutes per analysis, while the OMNIS NIR Analyzer delivers results in just a few seconds.

## EXPERIMENTAL EQUIPMENT

In this study, 17 tablets of paracetamol with varying water activity (0.23–0.85  $a_w$ ) were measured on an OMNIS NIR Analyzer (Figure 1) to create a prediction model for quantification. Samples were measured in reflection mode (1000–2250 nm) in 15 mm vials using a flexible holder and single-point measurement. The reference values were measured according to USP<922> Water Activity [3].



**Figure 1.** OMNIS NIR Analyzer Solid with 15 mm vial and Flexible holder OMNIS NIR.


**Table 1.** Hardware and software equipment overview.

Equipment	Article number
OMNIS NIR Analyzer Solid	2.1071.0010
Disposable vials, 15 mm, reflection	6.7402.110
Flexible holder OMNIS NIR	6.07402.300
OMNIS Stand-Alone license	6.06003.010
Quant Development software license	6.06008.002

## RESULT


The measured NIR spectra (Figure 2) were used to create a quantification prediction model for the water activity in paracetamol tablets. The quality of the prediction model was evaluated using the correlation diagram (Figure 3) which

displays a very high correlation between the NIR prediction and the reference values. The respective figures of merit (FOM) display the expected precision and confirm the feasibility during routine analysis.



**Figure 2.** Stacked NIR spectra of paracetamol tablets analyzed on an OMNIS NIR Analyzer Solid.

## RESULT WATER ACTIVITY IN PARACETAMOL TABLETS



**Figure 3.** Correlation diagram and the respective figures of merit for the prediction of water activity using an OMNIS NIR Analyzer Solid. The reference water content was determined using a Novasina LabMaster-aw neo according to USP<922>.

$R^2$	SEC ( $a_w$ )	SECV ( $a_w$ )
0.958	0.0278	0.0322

## CONCLUSION

This Application Note demonstrates the feasibility of determining  $a_w$  in paracetamol tablets quickly and easily. NIR spectroscopy offers users a fast, cost-effective, and highly

accurate alternative to other standard water activity measurement options. Additionally, NIRS analysis is non-destructive, completely reagent-free, and gives results in only a few seconds.

## REFERENCES

1. Pharmaceutical Trends: Water Activity Measurement - International Pharmaceutical Industry, 2021.
2. 〈922〉 Water Activity.  
[DOI:10.31003/USPNF\\_M12475\\_02\\_01](https://doi.org/10.31003/USPNF_M12475_02_01)
3. *USP 922 Water Activity Measurement - Novasina - Excellent new Method.*  
<https://www.novasina.ch/application/usp-922-water-activity/> (accessed 2024-08-27).

## CONTACT

Metrohm AG  
Ionenstrasse  
9100 Herisau

[info@metrohm.com](mailto:info@metrohm.com)

## CONFIGURATION



### OMNIS NIR Analyzer Solid 合固体和粘性品的近外光。

OMNIS NIR Analyzer 是一按照瑞士量准和生的近外光 (NIRS) 解决方案,用于整个生的常分析。使用新技和嵌入先 OMNIS Software 反在 NIR 光的速度、可操作性和活使用上。

OMNIS NIR Analyzer Solid 的点概:

- 在 10 秒以内量固体和粘性品
- 自化多位置量,即使在品不均匀,也能得可重的果
- 方便地嵌入自系,或者与其它分析技(滴定)
- 支持大量品容器



### 15 mm

123 个可封的玻璃一次性品瓶,直径 15 mm,用于分析反射中的固体。用于 XDS、DS2500 和 OMNIS 品系列的 NIR 固体分析。



### OMNIS NIR

可直径高 30 mm 的活支架,用于通反射品。

OMNIS  
A WHOLE NEW LEVEL OF PERFORMANCE

### OMNIS

允机版 OMNIS 件在一台 Windows™ 计算机上行。

特性:

- 可已含有一 OMNIS 可。
- 通万通可授平台行激活。
- 不可再外计算机上使用。

OMNIS  
A WHOLE NEW LEVEL OF PERFORMANCE

### Quant Development

用于在独立 OMNIS Software 安装套件中写和量化模型的件可。