

1 Purpose

This document describes the fixed bugs and problems as well as improvements in the following new software version in comparison with the previous version.

- New version: 1.2.55.1
- Old version: 1.2.55

2 New features

Topic	Description
Windows 11	Compatible with Windows 11. System requirements (in addition to the Windows 11 system requirements): <ul style="list-style-type: none">▪ Storage: 100 GB minimum, 500 GB recommended▪ USB connectors: USB 3.0 Note: Only use the provided Metrohm USB cable (6.021.08010), do not use third party USB cables. The usage of a powered USB Hub is recommended.▪ Mouse and keyboard: Required, USB, PS/2 style, or wireless▪ Screen resolution 1024 × 768 or higher▪ System backup: Network or local backup required for data archiving.

3 Improvements

Topic	Description
MCRL	Version 2 of the Metrohm Comprehensive Raman Library (MCRL) is added to the licensed library database.
Libraries	Supports the following library (since MIRA Cal DS 1.2.55): <ul style="list-style-type: none">▪ 66071634 MIRA Illicit & General Chemicals Library

4 Fixed bugs

None.

5 Known issues

- No support for Arabic language.
- Some language packs have errors or missing strings in the UI.
- One cannot view XTR processed data and baselined data for the same spectrum at the same time. One cannot view baselined data of an XTR processed spectrum.