

1 Purpose

This document describes features and known issues in MISA Cal software version 1.1.2.

2 Features

| Topic | Description |
|---|--|
| Available language packs | <ul style="list-style-type: none">EnglishChineseCzechFrenchGermanHungarianJapanesePolishTurkish |
| Spectral math functions | Weighted spectral addition and subtraction, Averaging, Standard deviation |
| Service support | Allow scripts to be run for service support: <ul style="list-style-type: none">Remove irregular pixelsRemove old licensed librariesTurn Bluetooth feature on / off |
| Viewing window data export | Modifications made to the viewing window (e.g., spectrum baseline, normalization, peak locations) can be exported in csv format |
| SERS Match | SERS match is shown in results screen |
| Custom Raman and SERS libraries creation | Ability to create, save and export custom Raman and SERS libraries |
| Custom Raman and SERS operating procedures creation | Ability to create, save and export custom Raman and SERS operating procedures |

