

1 Purpose

This document describes the fixed bugs in the following new software version in comparison with the previous version.

- New version: 4.2.300
- Old version: 4.2.200

2 New features

None.

3 Improvements

| Topic | Description |
|-------------------------------|--|
| Library support | MCUSP and MCRL libraries are now supported |
| Capitalize MIRA on login page | MIRA is now capitalized on the Login page |

4 Fixed bugs

| Topic | Description |
|--|---|
| IT Admin user has access to Samples page | The IT admin user no longer has access to the Samples page by simply exiting SST archive, Sync, or Service Portal pages |

5 Known issues

None.

6 Compliance

The current firmware version does not contain any modifications that affect conformity of Mira Cal regarding 21 CFR Part 11 and GAMP.