



772 Pump Unit – rinse

2.772.0130

Peristaltic pump for direct connection to Sample Processors. The "Rinse" version contains extensive accessories for rinsing.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 2.772.0130

Qt.	Order no.	Description
1 PCS	1.772.0110	772 Pump Unit
		Peristaltic pump for direct connection to Sample Processors.



1 PCS

6.1621.000

Container 10 L

As rinsing or waste container in automated systems.



2 PCS

6.1801.120

PVC tubing / 4 mm / 6 mm / 2 m



3 PCS

6.1805.060

FEP tubing / M6 / 60 cm

With light and kink protection



1 PCS

6.1808.170

Distributor for Sample Processors for rinsing and aspirating

Distributor for installation on Sample Processors for rinsing and aspirating when using externally connected pumps.



1 PCS

6.1812.000

PTFE tubing 4 mm / 6 mm / 4 m



2 PCS

6.1820.050

M6 / Olive adapter

For 6.1826.100 pump tubing, to peristaltic pumps.



1 PCS

6.1826.100

Pump tubing 6.4 mm / 1 m

Replaces the pump tubing used with the 772 Pump Unit and the 843 Pump Station (peristaltic).



1 PCS

6.1828.000

PVDF connection nipple

For 6.1621.000 canister.



3 PCS


6.2740.020

Spray nozzle

For the fine-spraying of the rinsing solution.



Optional accessories

Order no.	Description	
6.1826.160	Viton pump tubing For peristaltic pumps.	
6.2160.010	Adapter cable / 2 m / M8 socket / 2 x plug B Adapter cable / 2 m / M8 socket / 2 x plug B. For 772 Pump Unit.	