

843 Membrane Pump Station for Professional VA /CVS systems

2.843.0240

The 843 Membrane Pump Station enables the automatic rinsing and emptying of the measuring vessel in partly automated and fully automated Professional VA and CVS systems. It has two built-in membrane pumps that can be controlled via the remote signal.

The versions for VA and CVS are equipped with the complete accessories for the automatic rinsing and emptying of the measuring vessel.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 2.843.0240

|--|



The 843 Pump Station (membrane) has two built-in membrane pumps. These can be controlled directly via the interface using remote signals or manually by pressing a button.



4 PCS

6.1446.040

Threaded stopper / M6



1 PCS

6.1602.115

Bottle attachment / S 40 / green

For direct dosing from reagent bottles with S 40 thread (Merck...).



2 PCS

6.1621.000

Container 10 L

As rinsing or waste container in automated systems.





2 PCS

6.1819.010

PTFE tube

As gas inlet in washing bottle of VA Stand



1 PCS

6.1819.110

FEP aspiration tubing to canister

Aspiration tubing to canister attachment with M6 connector. For the production of eluents.



4 PCS

6.1820.010

Adapter for 1/8" gas line/M6 thread

For connecting a 1/8" gas line to an M6 thread.



For connecting up to 5 tubings



1 PCS

6.1829.070

6.1828.020

Tubing set for Professional VA/CVS instruments

Tubing set for the measuring head of Professional VA/CVS instruments. Comprised of 4 FEP tubings with M6 connector, tapering tip and protected against kink ($2 \times 82 \text{ mm} \log \text{und } 2 \times 110 \text{ mm} \log$).



1 PCS

6.2053.000

Cable clip

Cable clip for fastening cables and tubes



1 PCS

6.2621.000

Adjustable wrench

Max. opening: 20 mm. For IC instruments



Optional accessories

Order no.	Description
6.2141.200	Cable 732/819(remote)-733/820(fill/inject)
	Cable IC Detector (remote) - Advanced Separation Center (fill/inject)

6.2141.300 Remote cable 843 Pump Station to Sample Processor

Remote cable for the direct connection of the 843 Pump Station to the pump connectors of the Sample Processors



6.2148.010 Remote Box MSB

Additional remote interface for the connection of devices that can be controlled via remote lines. With permanently attached cable.

