



## Metrosep Carb 2 - 250/4.0

6.1090.430

The Metrosep Carb 2 - 250/4.0 IC column is particularly suitable for the determination of carbohydrates using alkaline eluents and pulsed amperometric detection. The high-capacity anion exchanger column is based on a styrene-divinylbenzene copolymer. It is stable in the range of pH = 0 - 14 and separates monosaccharides and disaccharides. It is also suitable for the analysis of sugar alcohols, anhydrous sugars, amino sugars, etc. The 250 mm version of the Metrosep Carb 2 separation column is optimized for complex separations.

Below, the accessories are grouped into Scope of delivery and Optional accessories.  
Please keep this printout at hand for ordering replacement material.  
These lists may be subject to change.

### Scope of delivery 6.1090.430

Qt.	Order no.	Description
-----	-----------	-------------

2 PCS





6.2744.060

Threaded stopper

For UNF 10/32. Stopper for IC, e.g. for the sealing of columns.



## Optional accessories

Order no.	Description	
6.1015.200	<p>Metrosep BO3 Trap 1 - 100/4.0</p> <p>Trap column for the removal of borate contaminants from the eluent. The Metrosep BO<sub>3</sub> Trap 1 - 100/4.0 is used mainly in carbohydrate analysis with hydroxide eluents.</p>	 A cylindrical chromatography trap column with a tan body and green end fittings. A label is visible on the side.
6.1015.300	<p>Metrosep CO3 Trap 1 - 100/4.0</p> <p>Trap column for the removal of carbonate traces in hydroxide eluents.</p>	 A cylindrical chromatography trap column with a tan body and green end fittings. A label is visible on the side.
6.1090.500	<p>Metrosep Carb 2 Guard/4.0</p> <p>The Metrosep Carb 2 Guard/4.0 effectively removes contamination, thus protecting the analytical separation column. The design of the guard column has been selected in such a way that its influence on the chromatographic separation can be ignored.</p>	 A guard column with a tan body and green end fittings. It has a label that reads 'Metrosep Carb 2 Guard', '4.0 µm', '100 Å', and '100 mm'. It is shown at an angle.
6.1090.510	<p>Metrosep Carb 2 S-Guard/4.0</p> <p>The Metrosep Carb 2 S-Guard/4.0 effectively removes contaminations, thus protecting the analytical separation column. The design of the guard column has been selected in such a way that its influence on the chromatographic separation can be ignored.</p>	 A guard column with a tan body and green end fittings. It has a label that reads 'Metrosep Carb 2 S-Guard', '4.0 µm', '100 Å', and '100 mm'. It is shown at an angle.