



# Reaction vessel long for stability measurements

6.1429.050

Reaction vessel for stability measurements involving volatile oils, mineral oils and biodiesel blends with the Biodiesel Rancimat. In combination with 6.2418.130. Set of 100 pieces.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

## Optional accessories

Order no.	Description
6.1808.090	Thread adapter M8 outer / M6 inner Outer thread M8, inner thread M6, e.g. for mounting 6.2418.XXX Air tubes in 6.1429.XXX Reaction vessels of stability measuring instruments.



6.2418.130

Air tube long for biodiesel measurements

Set of 100 pieces. In combination with 6.1429.050.

