



OMNIS Rod Stirrer Titration

2.1006.0020

Rod stirrer for connection to one OMNIS Titrator (cable length 1.1 m).

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.



These lists may be subject to change.

Scope of delivery 2.1006.0020

Qt.	Order no.	Description
1 PCS	1.1006.0020	Rod Stirrer Titration
		Rod stirrer for connection to one OMNIS Titrator (cable length 1.1 m).



Optional accessories

Order no.	Description	
6.01900.010	Stirring propeller 30 mm ETFE	
	Stirring propeller 30 mm ETFE with adjustable height for use with Rod Stirrer "Titrator" or Rod Stirrer "Sample Robot"	
<hr/>		
6.01900.030	Stirring propeller 20 mm ETFE	
	Stirring propeller, 20 mm, with adjustable height for use with Rod Stirrer "Titrator" or Rod Stirrer "Sample Robot"	
<hr/>		