



Double Pt sheet electrode

6.0309.100

Glass shaft electrode with two platinum sheets (0.15 x 8 x 8 mm) which are polarized for redox titrations (bivoltammetric titration). This electrode is well suited for

- Vitamin C determination
- Coulometric nitrogen determination
- Bromine index
- Sulfurous acid (SO₂) in wine
- Oxygen content according to Winkler

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.


Scope of delivery 6.0309.100

| Qt. | Order no. | Description |
|-----|-----------|-------------|
|-----|-----------|-------------|

Together with 6.2008.050 storage vessel holder, it provides a support for the electrode on 807 Dosing Units.



Optional accessories

| Order no. | Description |
|------------|---|
| 6.2104.020 | <p data-bbox="375 224 662 257">Electrode cable / 1 m / F</p> <p data-bbox="375 280 1157 358">For connecting electrodes with Metrohm plug-in head G to Metrohm instruments (socket F).</p>  |
