



# Consumable Kit, cylinder unit OMNIS special 20 mL

6.05700.451

Cylinder unit special 20 mL, complete with cylinder, piston, cylinder base, valve disk, centering tube, housing, and data chip. For 6.03004.220 OMNIS cylinder unit.

Below, the accessories are grouped into Scope of delivery and Optional accessories.  
Please keep this printout at hand for ordering replacement material.  
These lists may be subject to change.

## Scope of delivery 6.05700.451

Qt.	Order no.	Description
-----	-----------	-------------

0 PCS

6.2803.000

Grease (1 oz.)

Special quality (without silicone). For grease-free ground joint connections see 6.2713.XXX



---

1 PCS

6.2803.010

Grease (2 g)

Special quality (without silicone). For grease-free ground joint connections see 6.2713.XXX



## Optional accessories

Order no.	Description
6.03004.220	<p>Cylinder unit OMNIS special, 20 mL</p> <p>Intelligent 20 mL cylinder unit for one OMNIS Titrator, Titration Module or Dosing Module. This cylinder unit is especially recommended for the following solutions:</p> <ul style="list-style-type: none"><li>• Aqueous alkaline solutions</li><li>• Titrant 5</li><li>• Silver nitrate solutions</li><li>• Nonaqueous alkaline solutions</li><li>• Permanganate solutions</li><li>• EDTA solutions</li></ul> <p>Includes dosing tubing and antidiffusion tip.</p>

