



Metrosep C PCC 1 VHC/4.0

6.1010.320

Cation preconcentration column with very high capacity

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 6.1010.320

| Qt. | Order no. | Description |
|-------|------------|------------------|
| 2 PCS | 6.2744.060 | Threaded stopper |

For UNF 10/32. Stopper for IC, e.g. for the sealing of columns.

