



## Titration vessel lid with 5 openings

6.1414.010

Below, the accessories are grouped into Scope of delivery and Optional accessories.  
Please keep this printout at hand for ordering replacement material.  
These lists may be subject to change.

### Optional accessories

| Order no.  | Description  |
|------------|--|
| 6.1440.010 | Gas inlet and overflow tube with valve<br>For bubbling or covering of the sample by inert gas. |



6.2727.010

Intermediate SGJ sleeve

Used with reference and auxiliary electrodes in measuring and titration vessel with 6.1414.010 lid

