



807 Dosing Unit ETFE 5 mL

6.1575.150

807 Dosing Unit ETFE 5 mL for dosing alkali reagents or hydrofluoric acid (HF) with Dosinos. Less precision and reproducibility is to be expected in comparison with glass cylinders.

Below, the accessories are grouped into Scope of delivery and Optional accessories. Please keep this printout at hand for ordering replacement material. These lists may be subject to change.

Scope of delivery 6.1575.150

Qt.	Order no.	Description
-----	-----------	-------------

1 PCS

6.1446.040

Threaded stopper / M6



1 PCS

6.2803.010

Grease (2 g)

Special quality (without silicone). For grease-free ground joint connections see 6.2713.XXX



Optional accessories

Order no.	Description
-----------	-------------

6.1566.150 Cylinder element ETFE 5 mL

For dosing units used for dosing alkali reagents or hydrofluoric acid (HF) with Dosinos. Less precision and reproducibility is to be expected in comparison with glass cylinders.



6.1574.150 Cylinder element glass 5 mL

Glass cylinder element complete with cylinder, piston, cylinder base, valve disk, and centering tube. For 807 Dosing Unit 5 mL.

