



# Cannula for gas feeding

6.1819.080

For stability measuring instruments.


Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

## Optional accessories

Order no.	Description
6.0913.130	Measuring vessel cover for stability measuring instruments With built-in conductometric measuring cell.



6.1808.050

Tubing adapter olive / M8 outer

1 M8 outer thread and 1 tubing olive. E.g. for thermostat jacket of exchange units and stability measuring instruments.

