



2060 The NIR-R Analyzer

A472060010C

The 2060 *The* NIR-R Analyzer is the next generation of process spectroscopy instruments from Metrohm Process Analytics. With its unique and proven design from the inside out, it delivers accurate results every *10 seconds*. It can provide non-destructive analysis of liquids and solids directly in the process line or in a reaction vessel by using fiber optics and contact probes. It has been designed to connect up to five (5) probes and/or flow-cells. All five channels can be configured independently from each other using our versatile embedded proprietary software.

As part of the 2060 Platform, the 2060 *The* NIR-R Analyzer enables the unique separation of the human interface (HI) and NIR cabinet by optical fibers. This remote configuration allows the placement of both cabinets at different locations around the plant depending on customer preference and area classifications. Besides, it is available in three other versions: the 2060 *The* NIR Analyzer, 2060 *The* NIR-Ex Analyzer, and 2060 *The* NIR-REx Analyzer.

Below, the accessories are grouped into Scope of delivery and Optional accessories. Please keep this printout at hand for ordering replacement material. These lists may be subject to change.

Optional accessories

| Order no. | Description |
|-----------|-------------|
|-----------|-------------|

| | |
|-------------|----------------|
| 6.07410.000 | OMNIS NIR Lamp |
|-------------|----------------|

Spare lamp for 2060 The NIR Analyzer and OMNIS NIR Analyzer with built-in hour counter.

