



2060 The NIR-Ex Analyzer

A472060030C

The 2060 *The* NIR-Ex Analyzer is the next generation of process spectroscopy instruments from Metrohm Process Analytics. With its unique and proven design from the inside out, it delivers accurate results every *10 seconds*. It can provide non-destructive analysis of liquids and solids directly in the process line or in a reaction vessel by using fiber optics and contact probes. It has been designed to connect up to five (5) probes and/or flow-cells. All five channels can be configured independently from each other using our versatile embedded proprietary software.

Additionally, this analyzer is IECEx certified and fulfills the ATEX EU Directives. It has been designed with an approved purge/pressurization system together with intrinsically electronic devices, preventing any potentially explosive fumes or gases from the ambient air entering the analyzer enclosure. Furthermore, it is available in three other versions: the 2060 *The* NIR Analyzer, 2060 *The* NIR-R Analyzer, and 2060 *The* NIR-REx Analyzer.

Below, the accessories are grouped into Scope of delivery and Optional accessories. Please keep this printout at hand for ordering replacement material. These lists may be subject to change.

Optional accessories

Order no.	Description
-----------	-------------

6.07410.000	OMNIS NIR Lamp
-------------	----------------

Spare lamp for 2060 The NIR Analyzer and OMNIS NIR Analyzer with built-in hour counter.

